

# **Covid-19 Detection through X-Rays**

## **INTRODUCTION:**

In this pandemic which spread very quickly, it was difficult to test every person for the Covid-19 virus. Then came some variants which were not detectable early through RT PCR. Thus X-Rays and CT Scans were used to detect the virus. X-Rays were more popularly used to detect the Covid-19 virus.

As the cases rose at a very high rate during the 2nd wave, it was physically not possible to analyze each X-Ray and detect the Covid virus. Thus many Data Scientists and Radiologists came together to prepare different CNN models to detect Covid through X-Rays.

In this project, we also tried to create a model which will predict the Covid virus through X-Ray inputs.

## **Project Description:**

Deep learning has gained great importance in the field of Machine Learning. It has improved the performance of the model as well as increased the fields in which we can apply machine learning.

We have used Deep Learning to achieve our model. We have used Python inbuilt libraries in order to create our model.

Libraries used :

1. Keras
2. Matplotlib

3. Numpy
4. Pandas

The model is created is as follows:

Model: "sequential"

Layer (type)	Output Shape	Param #
conv2d (Conv2D)	(None, 222, 222, 32)	896
conv2d_1 (Conv2D)	(None, 220, 220, 64)	18496
max_pooling2d (MaxPooling2D)	(None, 110, 110, 64)	0
dropout (Dropout)	(None, 110, 110, 64)	0
conv2d_2 (Conv2D)	(None, 108, 108, 64)	36928
max_pooling2d_1 (MaxPooling2D)	(None, 54, 54, 64)	0
dropout_1 (Dropout)	(None, 54, 54, 64)	0
conv2d_3 (Conv2D)	(None, 52, 52, 128)	73856
max_pooling2d_2 (MaxPooling2D)	(None, 26, 26, 128)	0
dropout_2 (Dropout)	(None, 26, 26, 128)	0
conv2d_4 (Conv2D)	(None, 24, 24, 128)	147584
max_pooling2d_3 (MaxPooling2D)	(None, 12, 12, 128)	0
dropout_3 (Dropout)	(None, 12, 12, 128)	0
flatten (Flatten)	(None, 18432)	0
dense (Dense)	(None, 64)	1179712
dropout_4 (Dropout)	(None, 64)	0
dense_1 (Dense)	(None, 1)	65
Total params: 1,457,537		
Trainable params: 1,457,537		
Non-trainable params: 0		

The model needs to be trained on a dataset. Thus we used two sources to gain the dataset.


Source:

1. <https://github.com/muhammedtalo/COVID-19/tree/master/X-Ray%20Image%20DataSet>
2. <https://www.kaggle.com/fusicfenta/chest-xray-for-covid19-detection>

Now we used 178 images of COVID-19 affected X-Rays and 178 images of Normal X-Rays for training our model.


In order to validate our model we used 30 images of COVID-affected X-Rays and 30 images of Normal X-Rays.

The model trained for 100 epochs gave the following graph of training and validation accuracy:



After the training we validated our model and reached this confusion matrix. This confusion matrix shows that out of 30 Covid-19 affected X-Rays, 27

Covid-19 X-Rays were detected as Covid-19 affected X-Rays by our model while only 3 images were detected Normal.



Also we used an Image of Covid-19 affected X-Ray to predict whether the person has Covid-19 or not. And the model predicted it rightly as Covid-19 affected X-Ray.

```
In [16]: import numpy as np
# from google.colab.patches import cv2_imshow
import cv2
from keras.preprocessing import image
xtest_image = image.load_img('./Dataset/Prediction/ryct.2020200034.fig5-day7.jpeg', target_size = (224, 224,3))
xtest_image = image.img_to_array(xtest_image)
xtest_image = np.expand_dims(xtest_image, axis = 0)
results = model.predict_classes(xtest_image)
# training_set.class_indices
imggg = cv2.imread('./Dataset/Prediction/ryct.2020200034.fig5-day7.jpeg')
print("This Xray Image is of positive covid-19 patient")
imggg = np.array(imggg)
imggg = cv2.resize(imggg,(400,400))
plt.imshow(imggg)
# cv2_imshow(imggg)
# print(results)
if results[0][0] == 0:
    prediction = 'Positive For Covid-19'
else:
    prediction = 'Negative for Covid-19'
print("Prediction Of Our Model : ",prediction)
```

C:\Users\asus\anaconda3\lib\site-packages\keras\engine\sequential.py:450: UserWarning: `model.predict\_classes()` is deprecated and will be removed after 2021-01-01. Please use instead: `np.argmax(model.predict(x), axis=-1)`, if your model does multi-class classification (e.g. if it uses a 'softmax' last-layer activation).\* `(model.predict(x) > 0.5).astype("int32")`, if your model does binary classification (e.g. if it uses a 'sigmoid' last-layer activation).  
warnings.warn("`model.predict\_classes()` is deprecated and ")

This Xray Image is of positive covid-19 patient  
Prediction Of Our Model : Positive For Covid-19

## **Conclusion:**

We have created a model which can predict Covid-19 X-Ray with an accuracy  $> 90$ .